



Installing GASP in Medtech32 and Evolution

Note

This installation does not overwrite the existing installation – at the completion of the installation the practice will be able to access both the old and new Gasp form (until such time as the old form is inactivated).

Step 1 - Download the installation file

Go to <https://tools.comprehensivecare.co.nz/gasp-instructions-hhb> and click download for either Evolution or Medtech32

Step 2 - Install the new Advanced Form

Medtech32: Setup -> Advanced Forms -> Advanced Forms Manager

Advanced Form -> Import Disk

Evolution: File -> Options

Advanced Forms -> Advanced Form Manager

Actions -> Import Form

Step 3 – Testing

You should logoff and logon again before doing these tests.

Find a few patients that have an NHI, age, and ethnicity and verify that the content of the form looks reasonable and that there are no errors.